


THESE TWO CASES FIT ON THE PORT SIDE OF THE CABIN UNDER THE WINDOWS.

SECURE EACH CASE TO THE SIDEWALL WITH TWO STRAPS AND BUCKLES FROM THE KIT'S UNUSED CANVAS SEATS.


THE LETTERING AND STRAPS SHOW NICELY THROUGH THE WINDOWS.

## MEDICAL SUPPLY CASES FOR THE MODELART CH-34C MEDEVAC HELICOPTER

SMALL PACKS HERE



LARGE CASES HERE



THESE FOUR MEDICAL CORPSMAN'S PACKS FIT ON THE REAR PORT SIDE OF THE CABIN ABOVE THE WINDOWS.

LINE THEM UP ON THE FRAME RAIL (PART 1) AND SECURE THEM TO THE SIDEWALL WITH A LONG HORIZONTAL STRAP.

I PRINTED MY CH34 KIT SLIGHTLY UNDERSIZE (ABOUT 1/35 SCALE), ON LETTER-SIZED PAPER. THIS SHEET IS LAID OUT ON 8.5" X 11" PAPER. PRINT THIS SHEET SLIGHTLY OVERSIZE TO FIT A 1/32 SCALE MODEL.

FEEL FREE TO COPY & DISTRIBUTE THIS SHEET FREE OF CHARGE. COMMERCIAL REPRODUCTION PROHIBITED.

COPYRIGHT © 2003, DAVID SAKRISON, RIPON, WI USA  
ALL RIGHTS RESERVED.